

LCLS Single-Shot Relative Bunch Length Monitor System

-An Overview-

*M. Dunning
G. Travish
J. Rosenzweig
Particle Beam Physics Lab
UCLA Dept. of Physics
Los Angeles, CA 90095
October 11, 2005*

PB

mdunning@physics.ucla.edu

Intro: *The LCLS* Beamline




Linac Coherent Light Source

Stanford Linear Accelerator Center

Stanford Synchrotron Radiation Laboratory

<http://www-srl.slac.stanford.edu/lcls/>

- Operational 2009
- World's first X-ray FEL
- 1.6 cell S-band photoinjector
- 2 bunch compressors
- 100 m undulator



Intro: *The LCLS*

Relevant Parameters

	<u>1.0 nC</u>	<u>0.2 nC</u>	
Nominal electron energy, BC1	0.25	0.25	GeV
Nominal electron energy, BC2	4.3	4.3	GeV
Peak current	3400	2500	A
Nominal RMS bunch length, BC1	200	60	μm
Nominal RMS bunch length, BC2	20	8	μm
Nominal RMS bunch duration, BC1	667	200	fs
Nominal RMS bunch duration, BC2	67	27	fs
Max single bunch repetition rate	120	120	Hz

Intro: *The Problem*

- High-quality lasing: tight beam parameters
 - Longitudinal feedback systems needed
(along with other diagnostics and feedback systems)
 - Bunch length
 - Energy
- PBPL to build bunch length monitor system
 - System will consist of two grating polychromators, one at each bunch compressor
(explained later)

Intro: *Possible Solutions*

- Streak Camera
- Interferometer
- Electro-Optic Techniques
- RF Deflecting Cavity
- Polychromator (Spectrometer)
(more later)




Intro: *System Requirements*

- Only relative bunch length is needed- not absolute bunch length
- Need two bunch length monitors- one at each bunch compressor[1]
- Single-shot
- Non-invasive
- Maintenance free for several days
- Possibility to run at 120 Hz
- Single-shot measurement resolution: 1-2 % of nominal bunch length
- Long term signal drift: <2% over ~24 hours

[1] J. Wu et al., SLAC-PUB-11276, May 2005.

Intro: *Phase Feedback*



Observables

- Bunch length σ_z
- Energy E

Controllables

- Linac voltage V_{rf}
- Linac phase ϕ_{rf}

- LCLS longitudinal feedback: 2 bunch length loops
 - BC1 bunch length \rightarrow Linac 1 RF phase
 - BC2 bunch length \rightarrow Linac 2 RF phase

Possible Solutions

- **Streak Cameras**
 - + Single-shot
 - + Wide dynamic range
 - Limited by temporal resolution (~ 200 fs at best)
 - Trigger jitter




Hamamatsu "FESCA-200" (Femtosecond Streak Camera).
Temporal resolution: 200 fs.



Possible Solutions

- **Interferometers**
 - + Can be single-shot
 - + High temporal (frequency) resolution
 - + Compact
 - Narrow dynamic range
 - Complex




RadiaBeam Technologies BLIS (Bunch Length Interferometer System)

<http://www.radiabeam.com/products/diagnostics/blis.html>

Possible Solutions

- **Electro-Optic Methods**

- + Single-shot
- + Non-invasive (?)
- + Temporal resolution
- Not yet mature
- Require expensive femtosecond lasers



P. Bolton et al., SLAC-PUB-9529.
Transverse probe geometry produces a spatial image of the bunch.


Also see:

<http://www.rijnh.nl/users/berden/ebunch.html>



Possible Solutions

- RF Deflecting Cavities
 - + Single shot
 - + Femtosecond resolution
 - May require separate RF system
 - Invasive (destroy measured shot)




The UCLA 9-cell X-band standing wave deflecting cavity. Courtesy Joel England.



Possible Solutions

- Polychromators
 - + Single-shot
 - + Temporal resolution
 - + Robust
 - Require relatively expensive detector & vacuum system



Possible Solutions

Summary


	Single-shot	Non-Invasive	Good Temporal Resolution	Maintenance Free
Streak Camera	Y	Y	N	Y
Interferometer	Y	Y	Y	N
Electro-Optic	Y	Y (?)	Y	Y (?)
RF Deflector	Y	N	Y	N
Polychromator	Y	Y	Y	Y

Single-Shot Spectrometer

Bunch length monitor locations

- After 4th chicane magnet of BC1, BC2


BC1



Single-Shot Spectrometer Design

- Use CSR/CER from bunch compressor chicane magnets


- Vacuum port window
- Focusing/turning mirror
- Entrance slit
- Grating
- Off-axis parabola (line focus)
- Multichannel detector (linear array of cryogenically cooled bolometers)




Single-Shot Spectrometer

Bunch Distributions

BC1



BC2



- Smooth parabolic distribution

+ Simple CSR spectrum


- Wake-induced double-horn

- Complicated CSR spectrum



Single-Shot Spectrometer

Challenge: BC2 CSR Spectrum



CSR energy spectrum after BC2.
Black curve: double-horn distribution
Blue curve: Gaussian distribution
Red curve: step function

From J. Wu, et al., SLAC-PUB-11275, May 2005.

- Double-horn distribution complicates CSR spectrum
 - Similar to Gaussian below 4 THz
 - Stay below 4 THz

Single-Shot Spectrometer

Challenge: Detectors

BC1

- Frequency range: 150-500 GHz
- ~ 20 channels
- Easy, but big
 - large vacuum chamber
 - large optics
- InSb hot electron bolometers

BC2

- Frequency range: 1-4 THz
- ~ 20 channels
- More challenging than BC1
- Needs special filtering
- Thermal composite bolometers?
- Need to research more



Single-Shot Spectrometer


Challenge: Beamline Integration

- Low-loss vacuum port window over desired frequency range (Diamond?)
- Cryostats: liquid helium & nitrogen
 - Helium hold time (weeks?)
 - Closed-cycle nitrogen system (Sterling Engine?)
- Windowless enclosure for detector system

Single-Shot Spectrometer

BC1 Detector Assembly

- InSb hot-electron bolometers
- 10 liter cryostat
- Helium hold time: 4-6 weeks!




250 mm

20-channel linear array of InSb hot-electron bolometers, courtesy QMC Instruments.


Conclusion

Some work done so far...

Brookhaven CER work



UCLA built ATF compressor.



Simulated CSR spectrum from FieldEye, a post-processor of TREDI.



Brookhaven Si Bolometer for CER detection.[1]

Ref: G. Andonian, this workshop.

Conclusion

Workplan

- Simulate CR exiting vacuum ports of BC1, BC2 & arriving at detector
 - TREDI/FieldEye simulations
- Choose detector type
 - Finalize bolometer evaluations
 - SLAC to purchase
- Continue to study system
 - Windowless vacuum enclosure
 - Dynamic range (grating, *in situ* tuning)
 - Calibration methods
- Mechanical design & beamline integration with SLAC
 - CAD design work
 - Finalized by SLAC
- Test system (SPPS or APS Linac)